


bionet®

# Bioprocess lab and pilot equipment



*Enhancing Knowledge*



FO-BABY

**F1**

F2

F3

M1

M2

ROSITA & MARTA

bionet® F

I/O

Temp.

Digital  
Sensors

Level

Ext. Analog 1

Ext. Analog 2

Heating  
blanket

Agitator

Ext. Pump

Overlay

Sparger

Condenser  
return

Bioreactor  
return

Condenser  
supply

supply



# OUR NEW **F1** IS AN ADVANCED, RELIABLE AND VERSATILE MODEL



The F1 series include autoclavable bench scale stirred bioreactors up to 10L, an airlift bioreactor, and single-use stirred and rocking bioreactors with a patented 2D movement making the system unique in the market. The F1 is designed to meet the challenging and widely diverse R&D requirements and small-scale biomolecules production by using microbial and animal cells for biopharmaceutical-, food-, agricultural- and other biotechnological applications. Though a serially produced and standardized model, to combine the highest technological solutions for the common market demands, it is well thought for its expansion and customization towards a range of special requirements.

## **Elevate Your Bioprocess:**

Experience the **premium quality** of a globally recognized, standardized bioreactor, complete with swift worldwide delivery and expert service support.

## **Unrestricted from Your Initial**

**Decision:** Our modular Plug&Play design paves the way for seamless **expandability**, ensuring you're prepared for future enhancements and functionalities.

## **Effortless Operation and**

**Maintenance:** Discover an **ergonomically** designed bioreactor that incorporates cutting-edge technology and reliable innovations, making your daily tasks smoother than ever before.

## **Empower Your Process Control:**

Harness the potential of **advanced software** and instrumentation, providing you with a **user-friendly** interface to expand your knowledge and achieve precise process control.

## **Versatility Redefined with the F1:**

Tailor your bioprocessing needs with ease. Our high **configurability** accommodates various vessel volumes, bioreactor types, and combinations, offering unparalleled flexibility.

# CONFIGURABILITY

**Unrivalled product, groundbreaking technology:** The F1 unit provides the versatility to work with a wide range of applications and combine them, by maximizing the interchangeability of vessels. There are **several models** in Single or TWIN configuration, in multiple and single-use technologies, and in different **application versions**. **These include glass** stirred microbiology and cell culture, airlift, photobioreactor, and the CTB single-use **2D rocking patented** bioreactor. Each of them has unique features.




The CTB bioreactor features its distinctive patented 2D motion, which makes it also perfectly suited for bacteria and yeast fermentation processes.

Its expansion channel feature, also patented, allows a very wide working volume range, reducing the number of bioreactors in the chain and all the benefits this entails.



Equipped with a **gas module** with more than 21.000 combinations between number of mass flow controllers, bioreactor inlet, and gas flows to cover all needs.

**Freely assignable** pumps that allow to switch between operation modes and be linked to advanced controls from a unique configuration.



## EXPANDABILITY

**Perfect for R&D phases:** Ideal for those dynamic research and development phases where the requirements are constantly evolving and uncertainties abound. The Plug&Play modular design offers a strategic advantage by allowing you to address dynamic process requirements without the need for a substantial upfront investment



**The bBreath exhaust gas analyzer**, the module you need to carry out O<sub>2</sub>/CO<sub>2</sub> composition analysis of your bioreactor exhaust gas (i.e. metabolic indicators) and get the maximum of this data to enhance the quality and productivity of your process.

The **bScale** is the module you need for the connection of scales, from a range of scale precisions and brands used in lab and production settings.

**The Variable Speed Pump (bVSP)** is the module you need when expanding the number of freely assignable pumps to cover multiple additions or withdrawal/harvesting actions according to a number of feeding/bleeding profiles or advanced controls.




**On-Site expansion:** enjoy hassle-free modifications at your location. Our experts will handle the modification of actuators such as integrating new mass flow controllers in our expandable gas module, ensuring your equipment remains cutting-edge without the need for factory returns.

# ADVANCED AUTOMATION



The new F1 is powered by its inseparable and renewed friend **ROSITA 2.0**, Bionet's proprietary automation software for laboratory use, which allows for unparalleled automation flexibility and tight control over the processes and provides the user with ways to visualize, analyze and manage the data.

A sophisticated and comprehensive software solution that can rival the most advanced licensed products in the market while offering a more cost-effective alternative



## ROSITA 2.0

**Configurable control loops** between a number of sensors and actuators of choice, beyond the standard.

Programmable control loops via **equations**.

Configuration of **soft sensors** that can be visualized but also used for control.

**DASHBOARD**

Single 1

10:25:57 13/09/2023

**Bioreactor live**

pH 7,24

Temp 29,20 °C

Agit 0 rpm

DO 106,68 %

Foam

No batch running

**Soft sensors**

ss1 ID: 1 0 cells/ml

agracion ID: 2 100,00 rpm

+ Add Soft Sensor

**New bioreactor live for a precise visual representation on the current status of your process**

STIRRER

Basic control

**Personalize your dashboard to align with your specific work-flows, enhancing the efficiency of your process handling.**

**Controls**

pH		Temperature (°C)		Agitation (rpm)		DO (%)	
PV	7,24	PV	29,20	PV	0	PV	106,68
SP	5,00	SP	38,20	SP	1.200	SP	80,00

Gas 3_---		Pump 3 (ml/min)	
PV	0	PV	0
SP	0	SP	25,00

ALARMS CHECK: [BCU#1]

#1 T 0°C H 0% RAM 55% CPU 8% (1) BCU CONNECTED V2.0.2



# PREMIUM QUALITY



Our unwavering commitment to quality is reflected in our practice of incorporating top-tier instrumentation and components from the market, ensuring that our product achieves the highest standards of performance.



**Your go-to choice for small-scale or sample productions:** the F1 can be designed, constructed, and qualified in strict accordance with GMP guidelines, ensuring unparalleled quality and compliance.



# DATA SHEET

	AUTOCLAVABLE						
	STIRRED						
	MICROBIOLOGY				CELL CULTURE		
MODELS	F1-1 MB	F1-3 MB	F1-5 MB	F1-8 MB	F1-10 MB	F1-2 CC	F1-4 CC
VESSEL							
Total volume [L]	2.2	4.3	7.1	11.4	13.3	13.4	5.2
Working volume (max) [L]	1.3	3.0	4.8	8.0	10.0	1.85	4.0
Working volume (min) [L]	0.35	0.65 [1]	0.8	17 [2]	1.5	0.4	0.9
Vessel Material	Borosilicate glass						
Vessel frame & lid	SS 316						
GENERAL							
Total footprint on bench (HxWxD) [mm] - SINGLE	746 x 854 x 542		779 x 910 x 542		844 x 920 x 542	746 x 854 x 542	746 x 910 x 542
Total footprint on bench (HxWxD) [mm] - TWIN	746 x 1388 x 593		779 x 1500 x 593		844 x 1520 x 593	746 x 1388 x 593	746 x 1500 x 593
AGITATION							
Agitator	Top mounted - Single mechanical seal						
Impellers	2x Rushton						1x Pitched blade
Speed [rpm]	80-2000	80-1800	80-1500	80-1200	80-1000	80-500	
Motorpower [kW]	0.092	0.314				0.092	
GAS MODULE							
Gas lines (MFCs)	Standard: 1 / Optional: Up to 5						
Gases	Gas line 1: Air - N <sub>2</sub> / Gas lines 2 to 5: Air - CO <sub>2</sub>						
Gas inlet to vessel	Gas Line 1: Sparger / Gas line 2: Sparger / Gas lines 3 & 4: Selectable[4] - Sparger						
Gas flows	Low: 20-2000 sccm / Mid: 0.1 - 9 slpm / High: 0.1 - 9 slpm						
Condenser	•	•	•	•	•	•	•
High capacity condenser	-	-	•	•	•	-	-
Filter at exhaust gas	•	•	•	•	•	•	•
DOSAGE MODULE							
Pumps	Standard: 2 fixed speed pumps / Optional: 4 fixed speed pumps						
TEMPERATURE CONTROL							
Temperature range	10°C to 80°C						
Cooling	Circuit with automatic valves from external chiller to vessel jacket						
Heating	Electrical resistance						
INSTRUMENTATION							
Basic package	pH - Dissolved Oxygen - Temperature - Level						
EXPANSION POSSIBILITIES							
Scales (bScale)	•	•	•	•	•	•	•
Light Control Module (LCM)	-	•	•	-	-	-	-
Optical density	-	•	•	•	•	•	•
Exhaust Gas Analysis (bBreath)	•	•	•	•	•	•	•
Dissolved CO <sub>2</sub>	-	•	•	•	•	•	•
Viable Cells	-	•	•	•	•	•	•
UTILITY REQUIREMENTS							
CHILLED WATER							
Minimum flow rate	6 L/min SINGLE - 12 L/min TWIN						
Minimum pressure	0.8 barg						
COMPRESSED AIR SUPPLY	2-3 barg						
POWER SUPPLY							
	230 V (± 10 %), 50 Hz, Max power consumption 2500 W						
	120 V (± 10 %), 60 Hz, Max power consumption 2500 W						
	100 V (± 10 %), 50 Hz, Max power consumption 2500 W						
	Protection class IP 21						

[1] Standard is 1.5L minimum working volume. This 0.65L is with an optional Add-on Kit

[2] Standard is 2.65L minimum working volume. This 1.7L is with an optional Add-on Kit

[3] Depends on Draft Tube height


[4] Select from the software interface the bioreactor inlet for that gas line



			SINGLE-USE			
		AIRLIFT	STIRRED	ROCKING BIOREACTOR		
RE		CELL CULTURE				MICROBIOLOGY
C	F1-8 CC	F1-4 AL	F1-SU	CTB	CTB	
	9.65	4.7	3.0	6	10	
	7.0	[3]	2.4	2.5	5.0	
	1.0	[3]	1.0	0.2	0.2	
			Polycarbonate	Single-use bag		
			HDPE	-		
542 7	46 x 920 x 542	746 X 910 X 542	746 x 854 x 5427	46 x 1300 x 585		
x 593	746 x 1520 x 593	746 x 1500 x 593	746 x 1388 x 593	746 x 2220 x 593		
		Airlift effect	Top mounted	Rocking platform		
de		-	1x Pitched blade	-		
		-	80-500	0-35		
		-	0.092	-		
				Standard: 1 / Optional: Up to 4		
O <sub>2</sub> - CO <sub>2</sub> - N <sub>2</sub>				Gas line 1: Air - N2 / Gas lines 2 to 4: Air - O <sub>2</sub> - CO <sub>2</sub> - N <sub>2</sub>		
Sparger - Overlay / Gas line 5: Sparger - Overlay				Overlay		
: 0.2 - 18 slpm				25 - 2500 sccm		
	•	•	-	-	-	
	-	-	-	-	-	
	•	•	•	•	•	
nal: 2x integrated varibale speed pumps & 3x external variable speed pumps						
			Room temperature to 55°C	20°C to 40°C		
			-	Heat exchanger underneath the bag		
			Heating blanket	Heat exchanger underneath the bag		
el - Redox						
	•	•	•	•	•	
	-	-	-	-	-	
	•	•	•	-	-	
	•	•	•	-	-	
	•	•	•	-	-	
	•	•	•	-	-	
			-			
			-			
on 2500 W (UL compliant version) on 1100 W						


# AFTER-SALES SERVICE

**We acknowledge that acquiring our bioprocess equipment signifies the start of a lasting partnership, and rest assured, we'll be here to support you every step of the way!** Our core business is after-sales service, and this does not exist only in the form of maintenance and repair work but goes far beyond that. We offer product and process training to unblock gaps in processes or train your new employees.



**Our global after-sales service is a source of immense pride, delivered through our dedicated network of distributors who receive regular training at Bionet headquarters.**





Temp. Digital  
Sensors

Level

Ext. Analog 1

Ext. Analog 2

Heating  
blanket

Agitator

Ext. Pump

Overlay

Sparger

Condenser  
return

Bioreactor  
return

Condenser  
supply

supply

bionnet  
ENGINEERING  
32073



# Why Bionet?

A consolidated, configurable,  
scalable and integrated  
technological environment

---

A support team specialized in  
bioprocesses to help you from  
R&D to production

---

A partner for your entire project  
lifecycle, from initial technology  
configuration to after-sales



**bionet<sup>®</sup>**

**Bionet Engineering**

Parque tecnológico Fuente Álamo  
30320 Fuente Álamo, Murcia (Spain)  
Ph. +34 968 197 536  
Fax: +34 968 197 543

**[sales@bionet.com](mailto:sales@bionet.com)**  
**[www.bionet.com](http://www.bionet.com)**

**From  
Lab  
to  
Industrial**